




740SELECT monitor with Exergen TemporalScanner™ & cradle fully assembled

### 1 Remove the 740SELECT™ monitor and the accessory box from the shipping carton



1. Remove the plastic bag and monitor screen protector
2. Place the 740SELECT monitor and accessory box on a secure tabletop
3. The 740SELECT is shipped with the Li-ion battery installed

### 2 Open the 740SELECT standard accessory box and review contents

1. Identify the monitor part number located on the (monitor) rear label
2. Verify packout per table below

Part Number	NIBP	SpO <sub>2</sub>	Temp	ECG	Standard
01-02-1006	✓	✓	✓		✓
01-02-1007	✓	✓	✓		✓
01-02-1010	✓		✓		✓
01-02-1011		✓	✓		✓
01-02-1012		✓	✓		✓
01-02-1022	✓	✓	✓	✓	✓
01-02-1023	✓	✓	✓	✓	✓
01-02-1028		✓	✓	✓	✓
01-02-1029		✓	✓	✓	✓

#### Pack-Out Description

##### NIBP

BP start-up kit

##### SpO<sub>2</sub>

- SpO<sub>2</sub> patient cable
- SpO<sub>2</sub> sensor
- SpO<sub>2</sub> sensor start-up kit

##### Temp

• Exergen Temp cradle kit TAT-5000s-RS232-Qr Arterial, Deg C

- Quick reference guide addendum, Exergen

##### ECG

- ECG lead set
- ECG trunk cable

##### Standard

- Power supply
- Power cord
- User manual, addendums, and field service manual on CD
- Quick reference guide (monitor operation)

### 3 The following accessory item(s) may be placed in the monitor pack-out if purchased with the monitor

#### Masimo® CO<sub>2</sub>

- ISA™ CO<sub>2</sub> start-up kit
- IRMA™ CO<sub>2</sub> start-up kit

#### CO<sub>2</sub> Analyzer Accessory

- ISA CO<sub>2</sub> holster
- ISA CO<sub>2</sub> clamp adapter
- ISA CO<sub>2</sub> module holder

Note: The Microstream® Micropod™ start-up kit, cradle & clip (or cradle and VESA mount) are packaged separate from the monitor.

### 4 Patient connections


#### SpO<sub>2</sub>

The SpO<sub>2</sub> connector (left-hand side) and monitor input panel are keyed to assure proper orientation.

1. Align the SpO<sub>2</sub> connector and push into the input panel until the connection is seated.

- The connector is properly seated when you hear a click.

2. Complete the cable connection to the SpO<sub>2</sub> patient sensor.



Male NIBP Connector




#### NIBP Hose Connection CASMED MAXNIBP

1. Remove the BP hose from the BP start-up kit.
2. Align the female NIBP hose connector and push onto the male connector. You will hear a “click” when the connection is properly seated.
3. Gently pull on the BP hose to verify proper connection. Movement between the BP hose and connector is normal. The connection is designed to allow the hose to move slightly in and out and rotate while maintaining a secure connection.
4. Pull the metal ring to release the hose from the NIBP connector.

## 5 Installing the Exergen scanner cradle into the monitor


The cradle to hold the Exergen™ TempScanner must be attached to the Temp module.

1. Position the cradle so that the cradle tabs align with the Temp module.
2. Gently push the cradle fully into the Temp module slot.
  - The ears on the cradle tabs should engage and a click should be heard;
  - Gently pull on the cradle to ensure proper insertion of cradle into temp module.



## 6 Connect the (pre-attached) cable from the monitor to the Exergen TempScanner

1. Remove the bag with the scanner cable from the side of monitor, discard bag.
2. Align the cable and scanner.
3. Insert the cable into the bottom end of the scanner, push the cable until a click is heard.
4. Gently pull on the cable to ensure proper insertion of connector into the socket



## 7 Place the Exergen TemporalScanner into the cradle

1. Position the probe head of the scanner between the fingers of the cradle.
2. Push the scanner body into the cradle clips to secure.



Fully assembled cradle, cable & scanner secured in cradle

## 8 Connect the 740SELECT power supply and power cord

Power cord connected to power supply



Power cord Barrel connector

1. Connect the power cord to the power supply by fully inserting into the power supply.
2. Connect the barrel connector to the monitor.
3. Connect the power cord plug to a wall outlet.



## 9 Turn the 740SELECT Monitor on



1. Press the On/Off power button located on the front panel to turn the monitor on.
  - The On/Off switch illuminates when the monitor is turned on.
  - During boot-up a configuration setup test and POST is performed automatically (20-25 seconds).
  - Completion of the test and POST is confirmed by a two beep audio alert.
2. Verify the monitor battery charging LED is lit (when plugged into wall outlet).
  - Recommend fully charging battery before entering into service.
  - Reference 740SELECT user manual for operating instructions.