


# Nightingale MPC Technical Data

<b>General</b>	
Dimensions	8"W x 17.2"H x 17.5"D (203mm x 434mm x 445mm)
Power Requirements	90-135 / 180-265 VAC
Mains Frequency Range	47 – 63 Hz
Power Consumption	200W
Standards Conformance	UL 1950, EN 60950, EN 55022, EN 50082-1
<b>Device Markings</b>	
	European CE Mark according to Council Directive 93/42/EEC
<b>Environmental</b>	
Cooling	Fan
Operating Temperature	10° C to 35° C
Storage Temperature	-40° C to 65° C
Operating Humidity	>20% to <80% non-condensing
Storage Humidity	>20% to <80% non-condensing
Operating Altitude	0 to 3000 m
Storage Altitude	0 to 10600 m