



740SELECT™ CO₂ Accessories Ordering Guide: Covidien® Microstream® MicroPod™

The 740SELECT Series of multi-parameter patient monitors offers the option to continuously monitor End-Tidal Carbon Dioxide (etCO₂) using the Covidien Microstream MicroPod external etCO₂ module.

Covidien Microstream MicroPod




01-02-0841
Qty. 1

Covidien Microstream MicroPod External etCO₂ Module

• Plug and play Sidestream CO₂ analyzer. Microstream technology • Adult/Pediatric/Neonate

*Notes: Optional software required to enable CO₂ monitoring on the 740SELECT monitor;
Requires 740SELECT interface cable 01-02-0880*



01-02-0882
Qty. 1

Covidien Microstream MicroPod Cradle and Clip (RS09283)


Pole mounting option for MicroPod module

See individual items listed

01-02-0861
Qty. 1

Covidien Microstream MicroPod CO₂ Start-Up Kit


Kit contents: MicroPod etCO₂ module, FilterLine® Sample PAK and 740SELECT MicroPod interface cable



01-02-0883
Qty. 1

Covidien Microstream MicroPod Cradle and VESA Mount (RS09279)

VESA mounting option for MicroPod module



01-02-0880
Qty. 1

740SELECT Interface Cable (MicroPod)

Used to interface the Covidien Microstream MicroPod external etCO₂ module with a 740SELECT monitor

740SELECT CO₂ Software (Covidien Microstream)

The 740SELECT Open Technology Platform (OTP) allows CO₂ monitoring software to be installed at any time on any SELECT monitor as a CO₂ parameter upgrade.

01-02-0853
Qty. 1

Covidien Microstream MicroPod etCO₂ Software Factory Enable License

01-02-0856
Qty. 1

Covidien Microstream IPI Software Factory Enable License

Requires Covidien Microstream etCO₂ software factory enable license (PN 01-02-0853)

01-02-0855
Qty. 1

Covidien Microstream etCO₂ Software Field Enable License

01-02-0857
Qty. 1

Covidien Microstream IPI Software Field Enable License

Requires Covidien Microstream etCO₂ software field enable license (PN 01-02-0855)

Contact Zoe Medical customer support for ordering information at customersupport@zoemedical.com.

Product is  Marked

Dealer Imprint

740SELECT is a trademark of Zoe Medical. Microstream, Micropod, and FilterLine are trademarks of a Covidien company registered in the U.S. and foreign countries. Zoe Medical reserves the right to make changes to this sheet at any time without notice.