






## 740SELECT™ CO<sub>2</sub> Accessories Ordering Guide: Masimo® ISA™ & IRMA™ “Plug-in & Measure”

The 740SELECT Series of multi-parameter patient monitors offers the option to continuously monitor End-Tidal Carbon Dioxide (etCO<sub>2</sub>) using Masimo ISA (Sidestream) and IRMA (Mainstream) CO<sub>2</sub> analyzers.




### Masimo ISA CO<sub>2</sub> Sidestream Analyzer Monitoring Accessories

	01-03-0256 Qty. 1	<b>ISA CO<sub>2</sub> Sidestream Analyzer (800101)</b> • Plug-in and measure Sidestream CO <sub>2</sub> analyzer • Adult/Pediatric/Infant <i>Notes: Optional software required to enable CO<sub>2</sub> monitoring on the 740SELECT monitor; Requires 740SELECT interface cable 01-02-0494</i>
	01-02-0488 Box of 25	<b>ISA Nomoline Adapter (108220)</b> • Sampling line with female luer lock connector • With integrated water removal and hydrophobic bacteria filter • Adult/Pediatric/Infant • Multi-patient use
	01-02-0844 Box of 25	<b>ISA Nomoline Adapter with Male Luer Lock Connector (108210)</b> • With integrated water removal and hydrophobic bacteria filter • For general use • Length 2 meters • Single patient use
	01-02-0845 Box of 20	<b>ISA Nomoline Adapter Airway Adapter Set (108230)</b> • Sampling line with straight airway adapter • Adult/Pediatric • Single patient use • Integrated water removal and hydrophobic bacteria filter • 2 meters
	01-02-0848 Box of 25	<b>ISA T-Adapter (108250)</b> • Airway adapter with female luer lock connector • Adult/Pediatric • Single patient use
See individual items listed	01-02-0859 Qty. 1	<b>ISA CO<sub>2</sub> Start-Up Kit</b> <b>Kit contents:</b> ISA etCO <sub>2</sub> analyzer, Nomoline Airway adapter, Salter Labs CO <sub>2</sub> Nasal Cannulas (Adult, Pediatric, Infant, Intermediate Infant), 740SELECT interface cable
	01-02-0865 Box of 25	<b>Adult - Salter Labs® CO<sub>2</sub> Cannula Sample Line (4000)</b> • 7' (2.1 meters) supply tube, male luer connector • Single patient use • Latex free
	01-02-0867 Box of 25	<b>Pediatric - Salter Labs CO<sub>2</sub> Cannula Sample Line (4100)</b> • 7' (2.1 meters) supply tube, male luer connector • Single patient use • Latex free
	01-02-0869 Box of 25	<b>Infant - Salter Labs CO<sub>2</sub> Cannula Sample Line (4200)</b> • 7' (2.1 meters) supply tube, male luer connector • Single patient use • Latex free
	01-02-0871 Box of 25	<b>Intermediate Infant - Salter Labs CO<sub>2</sub> Cannula Sample Line (4300)</b> • 7' (2.1 meters) supply tube, male luer connector • Single patient use • Latex free
	01-02-0881 Qty. 1	<b>ISA CO<sub>2</sub> Analyzer Holster</b> • Attached to the left side of 740SELECT monitor • Provides a fixed or mobile solution for the ISA CO <sub>2</sub> analyzer


## Masimo ISA CO<sub>2</sub> Sidestream Analyzer Monitoring Accessories (continued)

	01-02-0935 Qty. 1	<b>ISA Analyzer Clamp Adapter</b> • Use with ISA analyzer module holder
	01-02-0936 Qty. 1	 <p><b>ISA Analyzer Module Holder</b> • Use with ISA analyzer clamp adapter • Movable jaw provides solution for rail attachment • ISA analyzer slides in and out of guide track</p>

## Masimo IRMA CO<sub>2</sub> Mainstream Analyzer Monitoring Accessories

	01-02-0843 Qty. 1	<b>IRMA CO<sub>2</sub> Mainstream Analyzer (200101)</b> • Plug-in and measure CO <sub>2</sub> analyzer • Adult/Pediatric/Infant <i>Notes: Optional software required to enable CO<sub>2</sub> monitoring on the 740SELECT monitor Requires 740SELECT interface cable 01-02-0494</i>
	01-02-0879 Box of 25	<b>Adult/Pediatric - IRMA CO<sub>2</sub> Airway Adapter (106220)</b> • <math>\leq 6\text{ ml}</math> deadspace • 0.3 cmH <sub>2</sub> O resistance @ 30 lpm • ET tubes > 4mm • Single patient use
	01-02-0877 Box of 10	<b>Infant - IRMA CO<sub>2</sub> Airway Adapter (106260)</b> • <math>\leq 1\text{ ml}</math> deadspace • 1.3 cmH <sub>2</sub> O resistance @ 30 lpm • ET tubes <math>\leq 4\text{ mm}</math> • Single patient use
See individual items listed	01-02-0860 Qty. 1	<b>IRMA CO<sub>2</sub> Start-Up Kit</b> <b>Kit contents:</b> IRMA etCO <sub>2</sub> analyzer, Adult/Pediatric Airway adapter, Infant Airway adapter, and 740SELECT interface cable

## ISA & IRMA Interface Cable

	01-02-0494 Qty. 1	<b>740SELECT Interface Cable (ISA/IRMA)</b> Used to interface an ISA or IRMA CO <sub>2</sub> analyzer with a 740SELECT monitor
---	----------------------	---

## 740SELECT CO<sub>2</sub> Software (Masimo)

The 740SELECT Open Technology Platform (OTP) allows CO<sub>2</sub> Monitoring Software to be installed at any time on any SELECT monitor as a CO<sub>2</sub> parameter upgrade.

01-02-0852 / Qty. 1

**Masimo etCO<sub>2</sub> Software Factory Enable License**

01-02-0854 / Qty. 1

**Masimo etCO<sub>2</sub> Software Field Enable License**

Contact Zoe Medical customer support for ordering information at [customersupport@zoemedical.com](mailto:customersupport@zoemedical.com).

Product is  Marked

Dealer Imprint

740SELECT is a trademark of Zoe Medical. ISA™, IRMA™, and Nomoline™ are Masimo, Inc. trademarks. Zoe Medical reserves the right to make changes to this sheet at any time without notice.