






740SELECT™ CO₂ Accessories Ordering Guide: Masimo® ISA™ & IRMA™ “Plug-in & Measure”

The 740SELECT Series of multi-parameter patient monitors offers the option to continuously monitor End-Tidal Carbon Dioxide (etCO₂) using Masimo ISA (Sidestream) and IRMA (Mainstream) CO₂ analyzers.




Masimo ISA CO₂ Sidestream Analyzer Monitoring Accessories

	01-03-0256 Qty. 1	ISA CO₂ Sidestream Analyzer (800101) • Plug-in and measure Sidestream CO ₂ analyzer • Adult/Pediatric/Infant <i>Notes: Optional software required to enable CO₂ monitoring on the 740SELECT monitor; Requires 740SELECT interface cable 01-02-0494</i>
	01-02-0488 Box of 25	ISA Nomoline Adapter (108220) • Sampling line with female luer lock connector • With integrated water removal and hydrophobic bacteria filter • Adult/Pediatric/Infant • Multi-patient use
	01-02-0844 Box of 25	ISA Nomoline Adapter with Male Luer Lock Connector (108210) • With integrated water removal and hydrophobic bacteria filter • For general use • Length 2 meters • Single patient use
	01-02-0845 Box of 20	ISA Nomoline Adapter Airway Adapter Set (108230) • Sampling line with straight airway adapter • Adult/Pediatric • Single patient use • Integrated water removal and hydrophobic bacteria filter • 2 meters
	01-02-0848 Box of 25	ISA T-Adapter (108250) • Airway adapter with female luer lock connector • Adult/Pediatric • Single patient use
See individual items listed	01-02-0859 Qty. 1	ISA CO₂ Start-Up Kit Kit contents: ISA etCO ₂ analyzer, Nomoline Airway adapter, Salter Labs CO ₂ Nasal Cannulas (Adult, Pediatric, Infant, Intermediate Infant), 740SELECT interface cable
	01-02-0865 Box of 25	Adult - Salter Labs® CO₂ Cannula Sample Line (4000) • 7' (2.1 meters) supply tube, male luer connector • Single patient use • Latex free
	01-02-0867 Box of 25	Pediatric - Salter Labs CO₂ Cannula Sample Line (4100) • 7' (2.1 meters) supply tube, male luer connector • Single patient use • Latex free
	01-02-0869 Box of 25	Infant - Salter Labs CO₂ Cannula Sample Line (4200) • 7' (2.1 meters) supply tube, male luer connector • Single patient use • Latex free
	01-02-0871 Box of 25	Intermediate Infant - Salter Labs CO₂ Cannula Sample Line (4300) • 7' (2.1 meters) supply tube, male luer connector • Single patient use • Latex free
	01-02-0881 Qty. 1	ISA CO₂ Analyzer Holster • Attached to the left side of 740SELECT monitor • Provides a fixed or mobile solution for the ISA CO ₂ analyzer


Masimo ISA CO₂ Sidestream Analyzer Monitoring Accessories (continued)

	01-02-0935 Qty. 1	ISA Analyzer Clamp Adapter • Use with ISA analyzer module holder
	01-02-0936 Qty. 1	 <p>ISA Analyzer Module Holder • Use with ISA analyzer clamp adapter • Movable jaw provides solution for rail attachment • ISA analyzer slides in and out of guide track</p>

Masimo IRMA CO₂ Mainstream Analyzer Monitoring Accessories

	01-02-0843 Qty. 1	IRMA CO₂ Mainstream Analyzer (200101) • Plug-in and measure CO ₂ analyzer • Adult/Pediatric/Infant <i>Notes: Optional software required to enable CO₂ monitoring on the 740SELECT monitor Requires 740SELECT interface cable 01-02-0494</i>
	01-02-0879 Box of 25	Adult/Pediatric - IRMA CO₂ Airway Adapter (106220) • $\leq 6\text{ ml}$ deadspace • 0.3 cmH ₂ O resistance @ 30 lpm • ET tubes > 4mm • Single patient use
	01-02-0877 Box of 10	Infant - IRMA CO₂ Airway Adapter (106260) • $\leq 1\text{ ml}$ deadspace • 1.3 cmH ₂ O resistance @ 30 lpm • ET tubes $\leq 4\text{ mm}$ • Single patient use
See individual items listed	01-02-0860 Qty. 1	IRMA CO₂ Start-Up Kit Kit contents: IRMA etCO ₂ analyzer, Adult/Pediatric Airway adapter, Infant Airway adapter, and 740SELECT interface cable

ISA & IRMA Interface Cable

	01-02-0494 Qty. 1	740SELECT Interface Cable (ISA/IRMA) Used to interface an ISA or IRMA CO ₂ analyzer with a 740SELECT monitor
---	----------------------	---

740SELECT CO₂ Software (Masimo)

The 740SELECT Open Technology Platform (OTP) allows CO₂ Monitoring Software to be installed at any time on any SELECT monitor as a CO₂ parameter upgrade.

01-02-0852 / Qty. 1

Masimo etCO₂ Software Factory Enable License

01-02-0854 / Qty. 1

Masimo etCO₂ Software Field Enable License

Contact Zoe Medical customer support for ordering information at customersupport@zoemedical.com.

Dealer Imprint

740SELECT is a trademark of Zoe Medical. ISA™, IRMA™, and Nomoline™ are Masimo, Inc. trademarks. Zoe Medical reserves the right to make changes to this sheet at any time without notice.