



Monitor Accessories

740SELECT™
Intelligent Monitoring Defined

SOLUTIONS ORDERING GUIDE

740SELECT Use Case Scenario

Accessory	Mobility Spot-Check Triage	Printer w/ Roll Stand	Mount to IV Pole	Transport Bed Rail Mounting	Wall Mounting	Ambulance Mounting	Fixed Surface	Barcode Reader w/ Roll Stand	Barcode Reader w/ IV Stand	Catalog #
Roll Stand	■									01-02-0172G
Roll Stand Accessory Kit	■									01-02-0821
Post Rail Attachment Kit			■	■						01-02-0819
12" Swivel Arm/Wall Channel					■					01-02-0265
VESA Adapter Plate					■	■	■			01-02-0809
Universal Swivel Mount						■	■			01-02-0173
Printer, 50mm		■								01-02-0807
Printer Mounting Kit		■								01-02-0820
Bar Code Reader								■	■	01-02-0815
C-Clamp and Holster									■	01-02-0816
C-Clamp, Holster & Flex Arm								■		01-02-0817



TRANSPORT

01-02-0948

740 SELECT Carry Case

- Use with all 740 SELECT Monitors



ROLL STAND

01-02-0172G

Roll Stand (5 wheels, 3-locking)

- Includes: Handle and wire basket with patient cable management hooks
- Requires Roll Stand Accessory Kit (CAS P/N 01-02-0821)

ROLL STAND ACCESSORY KIT



01-02-0821

740 SELECT Roll Stand Accessory Kit (required with CAS PN 01-02-0172G)

- Includes: Holster for power supply & power cable management (CAS P/N 01-02-0917), VESA Adapter Plate (CAS P/N 01-02-0809)

01-02-0917

Holster, Power Supply (attached to roll stand)

- Manages 740 SELECT power supply and power cord



01-02-0809

VESA Adapter Plate (includes 4 screws for attachment to monitor)

- Use with: Roll Stand (01-02-0172G), Swivel Arm (01-02-0265) or Universal Swivel Mount (01-02-0173)

POST RAIL ATTACHMENT KIT



01-02-0819

740 SELECT Post Rail Attachment Kit

- Kit Includes: Post Rail Clamp (01-02-0808) and L- Bracket (01-02-0813)
- Kit allows 740 SELECT to be attached to a bed rail or IV Pole

01-02-0808

Post Rail Clamp

- Use with L-Bracket (01-02-0813) to attach the 740 SELECT to a bed rail of IV Pole

01-02-0813

L- Bracket

- Use with Post Rail Clamp (01-02-0808) to attach the 740 SELECT to a bed rail or IV Pole



WALL MOUNTING

01-02-0265

12" Swivel Arm and Wall Channel

- Requires 740 SELECT VESA Adapter Plate (01-02-0809)

UNIVERSAL SWIVEL MOUNT



01-02-0173

Universal Swivel Mount

- Provides a fixed mounting solution for ambulances, table tops or a wall
- Requires mechanical attachment to a rigid structural surface
- 2.4" diameter base, 5" socket arm length
- Requires 740 SELECT VESA Adapter Plate (01-02-0809)

PRINTER



01-02-0807

External Printer, 50mm

01-02-0820

Printer Mounting Kit (Roll Stand)

- Kit Includes: Printer Mounting Bracket (01-02-0810) and Clamp Block (01-02-0918)

01-02-0810

Printer Mounting Bracket (Roll Stand)

- Requires Clamp Block (01-02-0918) to mount to roll stand



01-02-0918

Clamp Block

- Use with Printer Mounting Bracket (01-02-0810) to attach printer to the Roll Stand

01-02-0811

Printer Paper, Thermal, Pack of 12

POWER CORD & POWER SUPPLY



01-02-0396

U.S.A. Replacement Power Cord with Auto Lock and LED Plug

01-02-0806

Replacement Power Supply

MANUALS & GUIDES



01-02-1300

User Manual CD

21-22-0331

Field Service Manual (PDF)

21-22-0316

User Manual (PDF)

21-03-0337

Quick Reference Guide

BATTERY



01-02-0818

Rechargeable Li-Ion Battery Pack, 10.8V, 7800mAh

BARCODE READER



01-02-0815

Barcode Reader HS-1M (Medical Grade), with 14' coiled USB cable

01-02-0816

C-Clamp and Holster for IV Pole mounting of Barcode Reader (01-02-0815)

01-02-0817

C-Clamp and Holster with 12" Flex Arm, for mounting Barcode Reader (01-02-0815) to Roll Stand

BLOOD PRESSURE HOSE



01-01-0088

12' Coiled BP Hose with Female Bayonet Fittings

01-03-0255

6' Straight BP Hose with Female Bayonet Fittings

USB DRIVE



01-02-0814

Thumb Drive

INTERFACE CABLES

01-02-1100

Nurse Call Cable

01-02-0822

Serial Interface Cable, RJ45 to DB9 Female